

# Bioway Chemistry Reagent Series

## Cystatin C Test Kit


Detection of Cystatin C in Human Serum or Plasma on Chemistry Analyzers



Cat. No. R025K11

Cystatin C Test Kit

### SUMMARY OF TEST PROCEDURE



\*Refer to Figure 1 and the package insert for detail

**Table 1: Instrument Parameters\***

Calibration method	5 point non-linear	Slope of reaction	increase
Wavelength	Pr:340 nm Se: 700 nm	Sample volume	3 µl
Test method	2 point end	R1 volume	240 µl
Reaction temperature	37°C	R2 volume	80 µl

### INTENDED USE

**Bioway Chemistry Reagent Series Cystatin C Reagent Kit** (the Kit) is a latex-enhanced immunoturbidimetric assay intended for *in vitro* quantitative detection of Cystatin C in human serum or plasma on automated clinical chemistry analyzers.

### SUMMARY AND EXPLANATION

Cystatin C is a small cysteine proteinase inhibitor produced by all nucleated cells that can be freely filtered by the glomerular membrane and then nearly completely reabsorbed and degraded by the renal tubular cells. Thus, plasma concentration of cystatin C can be used as an indicator of glomerular filtration rate (GFR). Cystatin C levels are less dependent on age, sex, race and muscle mass compared to creatinine, making it a better assessment for kidney functions.

### TEST PRINCIPLES

The Kit utilizes latex-enhanced immunoturbidimetry to measure the Cystatin C level in human serum or plasma. During the test Cystatin C in the sample binds with the specific anti-Cystatin C antibody that is coated on latex particles to cause agglutination. The turbidity caused by agglutination is detected optically by chemistry analyzer. The change in absorbance is proportional to the level of Cystatin C in the sample. The actual concentration is obtained by comparing with a calibration curve with known concentrations.

### MATERIALS PROVIDED

#### Reagents:

<b>R1</b>	Glycine buffer solution
<b>R2</b>	Latex particles coated anti-Cys C antibodies

### MATERIALS NEEDED BUT NOT PROVIDED

- Automated chemistry analyzer.
- Cystatin C calibrator set and control set (available for purchase).

### INSTRUMENT

The Kit is applicable on most automated chemistry analyzers. Refer to specific instrument application for suggested settings.

### STORAGE AND STABILITY

Store the reagents at 2-8°C. Avoid direct sunlight. The Kit is stable through the expiration date when stored properly. R1 and R2 reagents are stable for 1 month at 2-8°C after opening.

### PRECAUTIONS

- The Kit is for *in vitro* diagnostic use only. Not for use in humans or animals.
- The instructions must be followed to obtain accurate results.
- Do not use the reagents beyond the expiration date.
- Treat all specimens as infectious. Proper handling and disposal procedures of specimens and test materials should be strictly followed.

### SPECIMEN COLLECTION AND HANDLING

Follow standard laboratory procedures to collect samples.

It is recommended to perform test immediately after sample collection. If the test cannot be done immediately, store sample at 2- 8° C for up to 2 days or at -20° C for up to 6 months. Avoid repeated freezing and thawing. **Do not use haemolysed samples.**

### TEST PROCEDURE (see Figure 1)

No pretreatment required for reagents and samples.

**Calibration:** 5 level calibrator set available for purchase. Recommend using Bioway calibrators for optimal results. Use multi-point non-linear calibration method.

**Test procedure:** see Figure 1 and Table 1 for instrument parameter setup. Refer to specific instrument application for suggested setting.

- Add 8 µl of sample and 240 µl of R1; mix well and incubate at 37°C for 180 seconds.
- Add 60 µl of R2, mix well and incubate at 37°C for 10 seconds.
- Take optical density measurement OD 1
- Incubate at 37 °C and take optical density measurement OD 2 at 480 seconds.
- Calculate  $\Delta OD = OD 2 - OD 1$

### RESULT

The Cystatin C value can be obtained by using the calculated  $\Delta OD$  to find the corresponding value on a calibration curve prepared with known values.

### EXPECTED VALUES

0.57 – 1.12 mg/L.

It is recommended for each laboratory to establish its own expected values.

### QUALITY CONTROL

Using Bioway controls with known concentration is recommended before each batch of tests to ensure the test is properly performed and all reagents and the instrument are functional as specified.

### LIMITATIONS

# Bioway Chemistry Reagent Series

## Cystatin C Test Kit

Detection of Cystatin C in Human Serum or Plasma on Chemistry Analyzers



1. The Kit is for *in vitro* use on automated chemistry analyzers only.
2. Hemolysis samples may cause inconsistent results.
3. The test result from the Kit should not be used as the only basis for definite diagnosis.

### PERFORMANCE CHARACTERISTICS

**Linearity:** 0 – 7.5 mg/L ( $R \geq 0.990$ )

**Precision:** Within Run:  $CV \leq 8\%$ ;  
Run-to-Run:  $CV \leq 10\%$

**Interference:** no interference detected for: Ascorbic acid (500g/L), Bilirubin (18 mg/dL), triglycerides (1000 mg/dL), hemoglobin (460 mg/dL), and RF (240 U/mL).

**Reagent Blank Absorbance:** at 546nm wavelength and 10 mm optical diameter, O.D.  $\leq 1.50$

### REFERENCES

1. A V Lewis, T J James, J B J McGuire and R P Taylor. Improved immunoturbidimetric assay for cystatin C. *Ann Clin Biochem* 2001; 38: 111 – 114.
2. Dharnidharka VR, Kwon C, Stevens G (August 2002). "Serum cystatin C is superior to serum creatinine as a marker of kidney function: a meta-analysis". *Am. J. Kidney Dis.* 40 (2): 221–6.
3. Kötting A, Selvin E, Stevens LA, Levey AS, Van Lente F, Coresh J (March 2008). "Serum cystatin C in the United States: the Third National Health and Nutrition Examination Survey (NHANES III)". *Am. J. Kidney Dis.* 51 (3): 385–94.
4. Mutsumi Tanaka, Kenje Matsuo, Masayasu Enomoto and Koji Mizuno. A Sol particle homogeneous immunoassay for measuring serum cystatin C. *Clin. Biochem.* 37 (2004) 27 – 35.
5. Roos JF, Doust J, Tett SE, Kirkpatrick CM (March 2007). "Diagnostic accuracy of cystatin C compared to serum creatinine for the estimation of renal dysfunction in adults and children--a meta-analysis". *Clin. Biochem.* 40 (5-6): 383–91.
6. Stevens LA, Coresh J, Schmid CH, *et al.* (March 2008). "Estimating GFR using serum cystatin C alone and in combination with serum creatinine: a pooled analysis of 3,418 individuals with CKD". *Am. J. Kidney Dis.* 51 (3): 395–406.

Not Intended for Sale in the United States.

Pointe Biotech (Nanjing) Co., Ltd.

No.85,Xingmin South Road, Jiang Ning District, Nanjing,P.R.China 211100

Tel:86-25-52425019,

Fax:86-25-52424836

[info@biowaydx.com](mailto:info@biowaydx.com)

