

# Bioway Chemistry Reagent Series UA Reagent Kit


Detection of Uric Acid in Human Serum, Plasma or Urine on Chemistry Analyzers



Cat. No. R055K11

UA Reagent Kit

## SUMMARY OF TEST PROCEDURE



**Table 1: Instrument Parameters\***

Calibration method	2-point Linear	Slope of reaction	Increase
Testing wavelength	Dλ : 520 nm Sλ : 600 nm	Sample volume	8 µl
Test method	1 point end	R1 volume	300 µl
Reaction temperature	37°C		

\*Refer to Figure 1 and the package insert for detail

## INTENDED USE



**Bioway Chemistry Reagent Series UA Reagent Kit** (the Kit) is an assay intended for *in vitro* quantitative detection of Uric Acid in human serum or plasma on automated clinical chemistry analyzers.

## SUMMARY AND EXPLANATION

Uric acid is a waste product derived from purine metabolism eliminating through the kidneys and gastrointestinal tract. High level serum uric acid may form monosodium urate crystals accumulated inside the bodies and causes gouty arthritis and renal diseases. It also indicates poor kidney and liver function, an increase in nucleoprotein metabolism, such as leukemia and polycythemia.

## TEST PRINCIPLES

The Kit utilizes enzymatic reactions to measure the amount of UA (µmol/L) in human serum or plasma.



Uric acid is oxidized in the presence of oxygen and uricase and the produced hydrogen peroxide then reacts with 4-aminoantipyrine (4-AAP) and N-ethyl-N-(hydroxy-3-sulfopropyl)-toluidine (TOOS) to yield quinoneimine which has the highest peak of absorbances at 520nm.

The decrease in absorbance at 520nm is directly proportional to the UA activity in the sample.

## MATERIALS PROVIDED

### Reagents:

R	Uricase Peroxidase 4-AAP TOOS Tris buffer, Ph7.7	300 U/L 1500 U/L 0.45 mmol/L 2 mmol/L 50 mmol/L
---	--	---

## MATERIALS NEEDED BUT NOT PROVIDED

- Automated chemistry analyzer
- UA control set (commercially available)

## INSTRUMENT

The Kit is applicable on most automated chemistry analyzers. Refer to specific instrument application for suggested settings.

## STORAGE AND STABILITY

Store the reagents at 2-8°C. Avoid direct sunlight. The Kit is stable through the expiration date when stored properly. The reagent is stable for 1 month at 2-8°C after opening.

## PRECAUTIONS

- The Kit is for *in vitro* diagnostic use only. Not for use in humans or animals.
- The instructions must be followed to obtain accurate results.
- Do not use the reagents beyond the expiration date.
- Treat all specimens as infectious. Proper handling and disposal procedures of specimens and test materials should be strictly followed.

## SPECIMEN COLLECTION AND HANDLING

Follow standard laboratory procedures to collect blood or serum preventing hemolysis. 10 times diluted urine is also acceptable. It is recommended to perform test immediately after sample collection. If the test cannot be done immediately, specimens may be stored at 2-8°C for 7 days.

## TEST PROCEDURE (see Figure 1)

Reagent 1 is liquid stable ready-to-use, no preparation needed.

**Calibration:** Recommend using included Bioway calibrator set for optimal results.

**Test procedure:** see Figure 1 and Table 1 for instrument parameter setup. Refer to specific instrument application for suggested setting.

- Add 8 µl of sample and 300 µl of R1; mix well and incubate at 37°C for 5 minutes.
- Take optical density measurement.
- Calculate average  $\Delta A$

## RESULT

The amount of UA in µmol/L can be obtained by the following calculation:

$$\text{UA } (\mu\text{mol/L}) = \frac{\Delta A_{\text{test}}}{\Delta A_{\text{standard}}} \times \text{standard solution } (\mu\text{mol/L})$$

Please refer to instrument application if testing under different conditions.

## EXPECTED VALUES

Male 208~428 µmol/L

Female 155~357 µmol/L

It is recommended for each laboratory to establish its own expected values

# Bioway Chemistry Reagent Series

## UA Reagent Kit

Detection of Uric Acid in Human Serum, Plasma or Urine on Chemistry Analyzers



### QUALITY CONTROL

Using commercially available controls with known concentration is recommended before each batch of tests to ensure the test is properly performed and all reagents and the instrument are functional as specified.

### LIMITATIONS

1. The Kit is for *in vitro* use on automated chemistry analyzers only.
2. The test result from the Kit should not be used as the only basis for definite diagnosis.
3. Samples with UA exceeding the maximum measurement range should be diluted with saline and retested.
4. Avoid high level Vitamin C contaminated in case of getting a lower test result.

### PERFORMANCE CHARACTERISTICS

**Linearity:** 0 – 1200  $\mu\text{mol/L}$  ( $R \geq 0.990$ )

**Accuracy:** Bias proportion 90%~110%

**Precision:** Within Run:  $CV \leq 4\%$ ;  
Run-to-Run:  $CV \leq 5\%$

**Interference:** no interference detected for: Unconjugated bilirubin ( $\leq 684 \mu\text{mol/L}$ ), Ascorbic acid ( $\leq 40 \text{mg/dl}$ ), Bilirubin ( $\leq 500 \text{mg/dl}$ ), Hemoglobin ( $\leq 3000 \text{NTU}$ ), Heparin, EDTA and Sodium Fluoride.

**Reagent Blank Absorbance:** At 520nm wavelength and 10mm optical diameter, O.D.  $\leq 0.10$ .

### REFERENCES

1. Fossati P. *et al.*, Clin. Chem., 26(2):227 (1980).
2. Morin L.G. *et al.*, J. Clin. Path., 60:691 (1973)
3. Praetorius E. *et al.*, J. Clin. Invest., 5:273 (1953)

Not Intended for Sale in the United States.

Pointe Biotech (Nanjing) Co., Ltd.

No.85,Xingmin South Road, Jiang Ning District, Nanjing,P.R.China 211100

Tel:86-25-52425019,

Fax:86-25-52424836

[info@biowaydx.com](mailto:info@biowaydx.com)

