

# Bioway Chemistry Control Series

## Cholyglycine Control Kit

1\*0.5 ml Liquid



Cat. No. R016C00

CG control Kit

### INTENDED USE

**Bioway Chemistry Control Series Cholyglycine Control Kit** is intended for the control of *in vitro* quantitative detection of serum Cholyglycine on automated clinical chemistry analyzers.

### PREPARATION

1. It is a ready-to-use control kit and no need for extra preparation before use.

### STORAGE AND STABILITY

Store the kit at 2-8°C till the expiry date. Avoid direct sunlight.

### PRECAUTIONS

1. The Kit is for *in vitro* diagnostic use only. Not for use in humans or animals.
2. The instructions must be followed to obtain accurate results.
3. Do not use the reagents beyond the expiration date.
4. In case of accidental ingestion, rinse the mouth and throat with excess water and immediately call for medical attention.
5. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local federal regulations.

### PERFORMANCE CHARACTERISTICS

Please see the values printed on the bottle label.

Not Intended for Sale in the United States.

Bioway LLC  
16776 Bernardo  
Center Dr. Suite 109  
San Diego, CA 92128  
Tel. (858) 487-0886  
Fax (858) 605-4208  
[info@biowayllc.com](mailto:info@biowayllc.com)

