

波音特生物科技（南京）有限公司

日立 (HITACHI 7170)

RBP

		Analyze	Calib.	Range	STD Conc	
Assay/Point	2Point End		10	16	34	0 0
Wave(Sub/Main)	700		340			
S. Vol. (Normal)	8.0		0.0	0		
S. Vol. (Dec.)	12.0		0.0	0		
S. Vol. (Inc.)			0.0	0		
Diluent	Water		0			
Reagent Vol. (R1)	225		0	#	0	
Reagent Vol. (R2)	0		0	00000	0	
Reagent Vol. (R3)	75		0	#	0	
Reagent Vol. (R4)	0		0	00000	0	
ABS. Limit	32000		Inc.	Twin Test		
Prozone limit	0	0	Lower			
Cell Det.	Detergent 1					

		Analyze	Calib.	Range	STD Conc
Calib. Type	SPLINE				
Point	6	Span Point			3
Weight	0				
Auto Calibration					
	Timeout			Chg. Over	
Blank	0	Change Lot			
Span	0	Change Bottle			
2Point	0				
Full	0				
SD Limit	999				
Duplicate Limit	500				
Sensitiv. Limit	0				
S1ABS Limit	-32000		32000		

		Analyze	Calib.	Range	STD Conc	
<Standard>		(1)	(2)	(3)	(4)	(5) (6)
Conc.		0.00	#	#	#	# #
Position		#	#	#	#	# #
Volume		8.0	8.0	8.0	8.0	8.0 8.0
<Pre-Diluent>						
Volume		0.0	0.0	0.0	0.0	0.0 0.0
Diluent		0	0	0	0	0 0
Calib. Code		0	0	0	0	0 0

“#” 由使用者自行输入浓度或位置。

地址：南京市江宁区科学园兴民南路85号
电话：025-52425019 传真：025-52424836-807